

RadPlanBio portal - User Documentation

27th June 2016

Contents

1 Overview of RadPlanBio

RadiationDosePlan- Image/Biomarker-Outcome-platform (RPB) is a collection of open source software systems, namely electronic data capture system (EDC - OpenClinica), medical image archive (PACS - Conquest) and patient identity management system (PIDG - Mainzelliste), which are integrated via RPB portal in order to deliver a core software infrastructure necessary to support translational radiotherapy research projects.




Figure 1: RPB portal home page.

1.1 Getting Started with RPB

1.1.1 Log In

To access the platform, in your web browser go to the RPB-portal home page. If you do not know the URL address, please contact the local RPB administrator at your site.


The image shows a login form titled "Login". Below the title, it says "RadPlanBio - portal login. Please provide your user credentials and choose a proper authentication service." There are three input fields: "Username", "Password", and "Authorisation method". The "Authorisation method" field is a dropdown menu with "OpenClinica" selected. Below these fields is a "Login" button, which is highlighted with a red rectangle.

Figure 2: RPB portal login page.

In order to navigate to the RPB login form, please click on first or second login link on the home page as it is highlighted in the Fig. 1. Afterwards you will be asked to provide your authentication credentials (user name and password) which will grant you role based access to the RPB features.

Note: RPB provides multiple authentication methods. The default way is to use the OpenClinica user account. With this account user can access complete range of RPB features (e.g. study as well as imaging data).

1.1.2 Change a Current Active Study

When you use RPB, you usually work within a specific Study and partner site. The name of the current active Study is above the main menu at left side of the screen. In order to change the active study you have to click on the name of the active study in top menu, see Fig. 3, which cause that the option **Change Study/site** will be shown. This option navigates to the study overview page, see :TP 2.1, where it is possible to change the active study for RPB (this change is also propagated to OpenClinica EDC tool).




Figure 3: Change user active study in RPB.

1.1.3 Read RSS news

RPB has a simple RSS module, where administrator can publish important news about development and maintenance of RPB platform. This news are shown at the RPB home page, to get there just click on the **Home** option in

main menu at the left side of the screen.

Note: there is also a direct link to the RSS source in the footer of RadPlaBio portal. This link can be used to subscribe for RSS news updated from standard RSS viewer software.

1.1.4 Log Out

To log out from RPB, click **Logout** (above the bread crumb navigation, on the right).

1.2 Portal Page Layout

The layout of RPB portal main page is separated to several areas with predefined features. Every page contains of the following areas depicted in Fig. 4: Top of page (Header) - a), Top menu (Below the header) - b), Main menu (Left side) - c), Body of page (d), Bottom of page (Footer) - e).




Figure 4: RPB portal page layout.

1.2.1 Top of Page

The top of the page is the header area of the page. It contains (from left to right):

- RPB partner site server
- Provider logo which navigates back to the RPB home page
- Notification of logged user with login/ logout and home button

1.2.2 Top Menu

Below the top page header there is a top menu. This menu consists of:

- Notification of current active study for user (left)
- Bread crumb navigation to show where in the hierarchy of the menu the user is currently (right)

1.2.3 Main Menu

The main menu is located at left side of the screen. It consists from the tree structured menu items which points to the concrete RPB modules. Which options are available depends on the type of logged user.

- Home

- CTMS - study management: clinical trial management module
 - Studies
 - Persons
 - Organisations
- EDC - data capture: electronic data capture module
 - Open EDC
 - Open Participate
 - Subjects/Events
 - Randomisation
- PACS - medical imaging: picture archiving and communication in medicine module
 - Subjects/Events/DICOM
 - DICOM Matrix
 - DICOM Lookup
- BIO - specimens: biobanking module
 - Subjects/Events/BIO
 - BIO Matrix
- PID - Patient identity: patient identity database module
 - Patient Search
- CMS - content management: content management system module
 - RSS Articles
- Data
 - Query
 - Import
- Administration: system administration module

1.2.4 Body of Page

What you see in the body of the page is depending on the module which you are accessing.

1.2.5 Bottom of Page

The bottom of the page contains:

- Home: Home page
- Platform: RPB platform web presence
- Collaboration: Map of partner sites defined this RPB portal installation
- Licence: List of RPB platform components with their licenses
- Software: List of software for download (RPB desktop client)
- Help: Link the on-line manual you are reading right now
- RSS: Direct link to RSS source of this RPB portal installation
- Impressum: Platform provider identification details
- Contact: Send an email to RPB portal administration
- Version Number: The version number is for the RPB portal software you are using

1.3 User Account

Each user of RPB platform is assigned a dedicated user account. Logged user is able to access his/ her user profile by clicking on user account name link in a header location of portal page. The profile view (Fig. 5) provides overview of user account details with assigned RPB user roles and filtered audit logs.

RadPlanBio - radiotherapy clinical research IT infrastructure

userdd

Logout

Home

client IP: 0:0:0:0:0:0:1
date: 27.06.2016

Home > User Account

User Account

User Roles

AuditLogs

?

Help

RPB Username:	userdd
OC Username:	userdd
Enabled:	<div>on</div>
Site:	Dresden
E-mail:	
LDAP:	<div>off</div>
API-key enabled:	<div>on</div>
API-key:	

Figure 5: RPB user profile page

2 EDC - data capture

EDC module provide the direct integration with OpenClinica EDC software used for the purpose of study data management. It shows and overview of studies the user is allow to access as well as study subject information. This module allows registration a new study subjects utilising patient identity management system to provide patient unique PID.

2.1 Studies Overview

In order to browse the studies in RPB system, the user has to navigate to the **Subjects/Events** menu item from the main menu of RPB portal, see Fig. 6.

The screenshot shows the RadPlanBio portal interface. At the top, the header includes 'Server: Dresden' and 'RadPlanBio - radiotherapy clinical research IT infrastructure'. The user is logged in as 'admindd' with a 'Logout' button. The client IP is '0:0:0:0:0:0:1' and the date is '27.06.2016'. The main navigation menu on the left includes 'Home', 'CTMS - study management', 'EDC - data capture', 'Open EDC', 'Open Participate', '★ Subjects/Events' (highlighted with a red box), 'Randomisation', 'PACS - medical imaging', 'BIO - specimens', 'PID - patient identity', 'CMS - content management', 'Data', and 'Administration'. The 'Active study/site: default-study' is shown at the top right. The main content area has tabs for 'Partner Sites', 'Subjects', 'Events', and 'Documents'. The 'Study' section shows a dropdown menu for 'Study' with 'Default Study' selected (highlighted with a red box). Below this is a table with columns: '#', 'Identifier', 'Site', and 'Commands'. The table contains one row with '# 1', 'Identifier default-study', 'Site Partner Sites', and 'Commands' (highlighted with a red box). The table has pagination controls showing '1' of '15' results. Below the table, it says 'There is one result' and 'Export' with a CSV icon. The footer includes links for 'Home', 'Platform', 'Collaboration', 'Licence', 'Software', 'Help', 'RSS', 'Impressum', and 'Contact', and the version 'radplanbio-portal version: 1.0.0.8 © 2013-2016 DKTK'.

Figure 6: Overview of studies available for logged in user in RPB.

The workflow for accessing studies via the RPB portal is as following:

1. Select the main **Study** from combo box (only studies which you have access to are shown).
2. Select a study **Partner Site** from list of the sites in data table.
3. (Optional) activation of the study, see below the sub section :TP 2.1.1.
4. Afterwards you can navigate to **Subjects** tab (sub section :TP 2.2).

2.1.1 Changing Current Active Study

In order to change the currently active study of the user. He has to use study overview feature and find a new study he would like to active (selecting main study definition from combo box and selecting the study site from data table). The Fig. 6 depict the button with the **Star** icon which is use to active a new study.

2.2 Study Subjects Overview

In this tab you will see an overview of **Subjects** which has been enrolled into the previously selected study and site, Fig. 7.

#	StudySubjectID	PID	Gender	Enrollment	Commands
1	xxa	DD-H22Y9VP2	m	2015-12-03	
2	xxb	DD-3K32A1MN	f	2015-12-03	
3	xxc	DD-Q235W02Y	f	2016-02-24	
4	xxx	xxx	m	2013-09-12	
5	xyy	0003Y0WZ	m	2013-09-13	
6	xxz	DD-0003Y0WZ	m	2013-09-18	

There are 6 results

Export

Figure 7: Study subjects.

2.2.1 Explanation of Possible ID Types

At this place it is important to explain that there are three types of patient identifiers which are used in RPB. Each of the type has its own purpose and one should understand what is the meaning of it:

- **StudySubjectID:** is identifying subject within one study. It can be automatically generated or manually specified (depending on study configuration). One patient taking part in multiple studies will have different StudySubjectID in each study. This ID has nothing to do with patient identity and cannot be use to obtain first or last name of patient no matter what.
- **PID** - pseudonym: is a unique identifier of the patient in RPB platform. One patient taking part in multiple studies will have the same PID in each one of them. It is generated from patient identity data (first name, last name, etc.) by patient identity management system in RPB. It is possible to not to use identity management system and manually provide PID, but in this case one has to guarantee the uniqueness of PID per patient and store the identity association somewhere outside of RPB. PID is also used in case of DICOM Patient ID flag when DICOM imaging data are uploaded.
- **SecondaryID:** sometimes it is useful to store patient ID from different systems and associate this ID to our study subject. And for this case we have a possibility to use this SecondaryID field.

2.2.2 Registering a New Study Subject

Before the new subject can be enrolled to the study, it is usually necessary to generate PID for this subject (When PID generator aka patient identity management system is used), see Fig. 8. The user is asked to fill some basic identity data. The optional identity data (like place or residence, ZIP code and birthname) are recommended if they are available for data entry personnel in order to guarantee uniqueness of generated PID. After data is provided the click on button **Generate** will trigger the process. The algorithm is trying to figure out if the patient with provided identity data is already in registered in the system (return his PID) or if it is new patient (generate new PID). After PID is provided the second **Submit** button is responsible for triggering the enrolment of the subject to the study.

Define a new Study Subject ✕

PID generator:	<input checked="" type="checkbox"/>
Study: *	Default Study
Partner Site: *	default-study
StudySubjectID: *	<input type="text"/>
Enrollment: *	<input type="text"/>
Gender: *	<input type="radio"/> Male <input type="radio"/> Female
Firstname: *	<input type="text"/>
Surname: *	<input type="text"/>
Birthname:	<input type="text"/>
Birth date: *	<input type="text"/>
Residence:	<input type="text"/>
ZIP:	<input type="text"/>
PID:	<input type="text"/>
<div>Generate</div>	
<div><div> Submit</div><div>✕ Reset</div></div>	

Edit selected Subject ✕

PID generator:	<input checked="" type="checkbox"/>
Study:	Default Study
Partner Site:	default-study
StudySubjectID:	xxb
SecondaryID:	
PID:	DD-3K32A1MN
Gender:	<input type="radio"/> Male <input checked="" type="radio"/> Female
Enrollment:	2015-12-03
	<div>Reidentify</div>
Firstname:	
Surname:	
Birthname:	
Birth date:	<input type="text"/>
Residence:	
ZIP:	

Register and enrol a new subject to the study.

Identify selected study subject.

Figure 8: Study subject commands.

2.2.3 Reading Study Subject Identity

Once the patient is successfully registered and enrolled into the study it is easy to get back his identity details when necessary (according to PID), see Fig. 8. The identity information stored in patient identity management system is retrieved when the user click on **Reidentify** button.

2.3 Scheduled Study Subject Events Overview

Tab **Events** is showing the table of events which have been scheduled for selected subject in this study. No more functionality is implemented here right now.

2.4 Direct Access to OpenClinica

By clicking on **Open EDC** link in the main menu at the left side of the page, you can directly navigate to OpenClinica EDC system.

2.5 Direct Access to Participate

By clicking on **Open Participate** link in the main menu at the left side of the page, you can directly navigate to RPB participate system for mobile data entry.

2.6 Randomisation

RPB portal randomisation module can be used to randomise subjects into treatment arms according to defined set of stratification criteria. The module is available only for those clinical trials which have been configured as randomised.

Server: Dresden RadPlanBio - radiotherapy clinical research IT infrastructure

adminidd Logout Home

client IP: 0:0:0:0:0:1
date: 27.06.2016

Active study/site: default-study

Menu

- Home
- CTMS - study management
- EDC - data capture
 - Open EDC
 - Open Participate
 - Subjects/Events
 - Randomisation**
- PACS - medical imaging
- BIO - specimens
- PID - patient identity
- CMS - content management
- Data
- Administration

Subjects

Reload Print Help

#	StudySubjectID	PID	Arm	Commands
1	xxa	DD-H22Y9VP2	Treatment Arms	>>
2	xxb	DD-3K32A1MN		>>
3	xxc	DD-Q235W02Y		>>
4	xxx	xxx		>>
5	xxy	0003Y0WZ		>>
6	xxz	DD-0003Y0WZ		>>

There are 6 results

Export

CSV

Figure 9: Randomisation module overview.

2.6.1 Randomise Study Subject

Once you enter the randomisation module you will be able to see the list of currently registered study subject within your active study as displayed in Figure 9. The user can distinguish whether the subject was already randomised by examining his treatment arm property.

For non randomised subjects the command button at the right side of each subject data row will trigger dialogue allowing to enter necessary stratification criteria for the randomisation process, see Figure 10.

After clicking on assign random treatment arm button, the randomisation takes place and assign one of possible study treatment arms to the subject. After this step the treatment arm for subject is stored within RPB. In the current version of RPB it is however still necessary to manually copy this information into OpenClinica EDC system. Within OpenClinica there are two solutions for how the treatment arm of the patient can be stored:

1. Treatment arm as subject group class: OpenClinica EDC system provides the ability to defined subject group classes. During setup of a concrete study one can declare treatment group class with treatment arm options and use this data element to store the randomisation result, see Figure 11.
2. Treatment arm as eCRF item field: Randomisation result can be also stored within one of the eCRFs of a specific study. During a setup of the study one has to create eCRF field element which can store the randomised treatment arm option, see Figure 11.

The choice which option to use (even both of them could be used together) depends on a concrete setup of randomised clinical trial.



Randomise Subject ✕


StudySubjectID:	xxa
Partner Site:	Dresden
Stratify Partner Site:	true
Arm:	

T: * Select One ▼
strata

Assign random treatment arm

Figure 10: Randomisation dialogue with stratification details.

Assign Subject to Group



Study Subject ID: * indicates required field.

Secondary ID:


Enrollment Date: 📅

Subject Group Class: Treatment group -- ▼

Notes:

Confirm Changes
Cancel

Randomised treatment arm as study subject group in EDC.



⏪
I Ein- ... (0/19)
II Anga... (0/20)
-- Select to Jump -- ▼
⏩

Titel: Ein- und Ausschlusskriterien

Instructions: Bitte alle Felder ausfüllen!

Page: Save Exit 📄

Gruppe

Behandlung Gruppe: (Wählen Sie eine Option) ▼
📄

Randomised treatment arm as EDC eCRF field.

Figure 11: Treatment arm in OpenClinica EDC.

3 PACS - medical imaging

RPB has PACS module which provides integration with research PACS system where all the study related DICOM data are uploaded with RPB desktop upload tool. RPB uses following information linking in order to store association between imaging data/ clinical study data/ and patient identity:

- DICOM patient ID: during upload of DICOM data the pseudonymisation step is performed which also ensures that all DICOM patient ID tags in original file are replaced with selected RPB patient pseudonym.
- DICOM study UID: the old study UID are replaced with randomly generated new study UID.

These two types of identifier are saved in CRF for specific study event with DICOM eCRF to ensure that we can always find appropriate study related imaging RPB research PACS.

Note: this replace and storing of DICOM ID is done automatically via the RPB desktop upload tool.

Once imaging data is successfully stored within RPB the user can navigate them by invoking **DICOM studies** menu item from main menu. Here the user can browse DICOM studies, see Fig. 12 or study series, of selected patient in currently active study. This DICOM data are able for download or viewing in integrated DICOM viewer.

Server: Dresden RadPlanBio - radiotherapy clinical research IT infrastructure

adminidd Logout Home

client IP: 10.44.89.97
date: 27.06.2016

Active study/site: default-study

Menu: Home, CTMS - study management, EDC - data capture, PACS - medical imaging, **Subjects/Events/DICOM**, BIO - specimens, PID - patient identity, CMS - content management, Data, Administration

Subjects Events **DICOM**

Drucken Hilfe

xxxy-[0003Y0WZ] /ST Imaging[1]

1 15 Spalten

#	UID	Description	Modality	Time	Kommandos
1	2.25.1922	helical	CT	000000	View Download
2	2.25.4763		RTSTRUCT		View Download
3	2.25.7527		RTPLAN		View Download
4	2.25.7652		RTDOSE		View Download

Es gibt 4 Ergebnisse

2 15 Spalten

#	UID	Description	Modality	Time	Kommandos
1	2.25.1027	1.25mm	CT	124506	View Download
2	2.25.1349	1.25mm	CT	124424	View Download
3	2.25.3774	Scout	CT	124325	View Download
4	2.25.4561	Injection Report	CT	124325	View Download
5	2.25.5387	Dose Report	CT	124325	View Download
6	2.25.7085	2.5mm	CT	124506	View Download
7	2.25.7746	2.5mm	CT	124424	View Download

Es gibt 7 Ergebnisse

Figure 12: Overview of DICOM studies for subject within study and event.

4 PID - patient identity

RPB PID module is used to safely access the information from RPB patient identity management system. In section :TP 2.2.2, it was described how a new PID pseudonym is generated and patient is enrolled as a subject of specific clinical trial.

4.1 PID - explanation

PID (patient identifier, person identifier, pseudonym) is a unique identification string which is used across the sub-systems in RPB in order to pseudonymise the patient identity. This means that patient identity is stored only in one place, so called, patient identity database. On the other hand when we want to associate some data in RPB to the specific patient (it can be clinical study data from CRF or DICOM data) we use the PID which uniquely represent the patient instead of real identity data.

4.2 Getting PID of registered patient

From time to time the data entry personnel knows the identity of the registered patient and need obtain PID pseudonym associated to this patient while it was registered. This is possible by navigation to the **Search patient PID** menu item in the main menu as it is displayed in Fig. 13.

Server: Dresden RadPlanBio - radiotherapy clinical research IT infrastructure

admindd Logout Home

client IP: 0:0:0:0:0:0:1
date: 27.06.2016

Active study/site: default-study

Menu

- Home
- CTMS - study management
- EDC - data capture
- PACS - medical imaging
- BIO - specimens
- PID - patient identity
 - Patient Search**
- CMS - content management
- Data
- Administration

Patients Studies DICOM

Search mode

IDAT: on

Search criteria

Firstname: *
Surname: *
Birth date: *
Residence:
ZIP:
Search Reset

New Reload Reidentify Print Help

#	PID	Surname	Firstname	Birth date	Commands
No matching result					

Figure 13: Lookup of PID when we know the identity information of registered patient.

5 How To

This section provides some generic but useful information about a specific usage of RadPlanBio.

5.1 DICOM upload for dummy patient (step-by-step guide)

Some studies in RadPlanBio require for the partner site a testing DICOM imaging or treatment plan upload before it can join the real study data collection. For such reason every partner site has an access to the **Demo Study** where such testing procedures can be performed. The whole process of DICOM upload for a dummy patient consists of three steps which are described in details below:

- Patient registration and enrolment into a study
 1. Login with your user credentials into RadPlanBio portal.
 2. On the left side (tree menu) click on **EDC - clinical study – Studies/Subject/Events** item to access an overview of studies that are available for you.
 3. From **Study** (combobox) select the study with name **Demo Study** (this is where the dummy patient will be created).
 4. From the data table which will be shown below select your partner site by clicking on appropriate data table row.
 5. Click on **Study Subjects** (tab) to get an overview of subjects already enrolled into a study.
 6. Click on the **New** (button) to register a new subject.
 7. Enter all the required dummy subject fields:
 - Study subject ID: you can use something like [PartnerSiteIdentifier]-DUMMY-001 (e.g. DD-DUMMY-001).
 - Enrolment date: e.g. today.
 - Gender: select one.
 - First name: dummy.
 - Last name: patient.
 - Birth date: choose any.
 8. Click on the **Generate** button to obtain unique patient pseudonym – it will be shown in the generated patient PID field.
 9. Click on the **Submit** button to enrol pseudonymised patient into a study.
- Scheduling the imaging study event
 1. From the left menu in RadPlanBio portal website select **EDC - clinical study – Open data capture tool**.
 2. EDC system (OpenClinica) website will open.
 3. In **Subject Matrix** you should be able to see subject enrolled into a currently active study. It should be Demo Study: [PartnerSite]. If it is not you should change the active study by clicking on the **Change Study/Site** link in upper menu.
 4. Click on empty document icon for specific subject row in DE Imaging column. And click on schedule to schedule this study event (data collection is only possible for scheduled events).
 5. Fill up some details for event (date) and schedule event for this subject.
- Uploading treatment plan (or different DICOM imaging) for registered subject with scheduled imaging study event:
 1. Start the RadPlanBio client application and login with your user credentials.
 2. Click on the button **Upload DICOM data....**
 3. From the **Study** (combo box) select the **Demo Study**.
 4. From data table displayed below select your site.
 5. Go to the next **Study Subject** (tab).
 6. Select the **dummy** subject you have created in the first step.
 7. Go to next **Study Events** (tab).

8. Select the scheduled **DE Imaging** study event.
9. Go to the next **DICOM** (tab), here select the element **Treatment plan** (or if you need to upload another DICOM modality select the appropriate element).
10. Go to the next **Summary** tab.
11. Click on the button **Upload DICOM** data.
12. Navigate to the folder on your computer where the treatment plan (or another DICOM study) you want to upload is stored. Select the folder and click **OK**.
13. Files within folder will be analysed and the structure of underlying DICOM study will be displayed in a tree view. Select the whole study or all necessary DICOM modalities (for treatment plan upload at least CT, RTDOSE, RTPLAN, RTSTRUCT modalities are required. When selection is done click **OK**.
14. Treatment plan is analysed and checked for consistency. The next dialogue will show you how treatment plan will be de-identified (patient pseudonym generated in the first step will be used).
15. You may want to click on buttons to preserve the descriptions of DICOM study and series. Otherwise these descriptions are going to be removed during the de-identification process.
16. When this is done click **OK**, the question will be shown to ensure that study and series descriptions do not contain any patient identity data.
17. Next you will be asked to harmonise the naming of organs contours from provided RTSTRUCT series.
18. Each named contour should be mapped with appropriate standard name. For specific contours (e.g. GTV, PTV, ...) it is possible to provide an additional information. Dual organs it is possible to choose between left/right in additional info combo box.
19. After original contours names are mapped to standard names click **OK**.
20. Question is displayed to ensure about the correctness of provided mapping, after click on **yes** the pseudonymisation followed by the upload of data procedure will start.
21. After upload is finished the **Data transfer was successful** message will be displayed.